

BE IN THE MOMENT

WITH HANSATON FOKUS
HEARING AIDS



A Sonova brand

 **HANSATON**
hearing & emotions

RECHARGEABLE TECHNOLOGY:

Energy for the whole day

Four HANSATON FOKUS AQ hearing aid styles with advanced Lithium-ion technology allowing

- > Up to 24 hours runtime between charges
- > 3 hours charging time to a full charge
- > To answer the phone with a simple tap on the ear (except for Performance Level 1)

24

Rechargeable Runtime (charging from 0 % to 100 % in less than 3 hours)

Up to 24 hours: Without streaming

Up to 16 hours: AirStream technology (TV Connector/PartnerMic)

Up to 11 hours: Bluetooth® phone and music streaming

HANSATON FOKUS

HANSATON FOKUS hearing aids come with a modern, unique design and the exceptional FOKUS sound performance. They're available in five performance levels powered by the innovative HearIntelligence™ technology. They also offer connectivity to all mobile phones¹ and hands-free phone calls.



Sterling Silver

Metallic Anthracite

Carbon Black

Cashmere Beige²

Mocha Brown²



¹ Smartphones and traditional mobile phones with compatible Bluetooth® profiles.
² Except AQ sound FS S

WIRELESS CONNECTIVITY:

Comfort for every day

All HANSATON FOKUS hearing aid styles provide:



Hands-free phone calls with all mobile phones¹




Enable your clients to stream videos and audio from their mobile devices directly to their hearing aids.



Explore optional accessories: TV Connector, PartnerMic, Remote Control, and Roger™ microphones.



Personalize the listening experience with the stream remote App



Take a tour of the stream remote App

¹ Smartphones and traditional mobile phones with compatible Bluetooth® profiles.
The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc.




HANSATON FOKUS: RIC hearing aids



HANSATON FOKUS: BTE hearing aids

“The great sound performance and connectivity you want, in a design you will love.”
CLIENT BENEFIT



- sound FS R312**
- Smallest HANSATON hearing aid
 - Powered by a 312 replaceable battery



- AQ sound FS R**
- HANSATON's smallest rechargeable hearing wonder
 - **Compatible with:** Comfort Charger Combi and HANSATON Charger



- AQ sound FS RT**
- Telecoil-equipped
 - **Compatible with:** HANSATON T Charger



- AQ sound FS S**
- Slim and stylish
 - **Compatible with:** HANSATON S Charger



- beat FS RT675 UP**
- HANSATON's most powerful BTE
 - Dual receiver with high output and gain
 - Telecoil-equipped
 - Powered by a 675 replaceable battery




- beat FS RT13 SP**
- Super Power BTE
 - Telecoil-equipped
 - Powered by a 13 replaceable battery



- AQ beat FS R**
- Dual receiver with high output and gain
 - **Compatible with:** Easy Line Charger Combi BTE and HANSATON Charger



- beat FS RT312 M**
- Moderate Power BTE
 - Telecoil-equipped
 - Powered by a 312 replaceable battery




— Open Dome / Cap Dome
 ... Vented Dome
 --- Power Dome / SlimTip


Receiver type	S	M	P	UP ³
Output / Gain (2cc)	111 / 47	114 / 51	122 / 59	132 / 71

³ Not compatible with AQ sound FS S


beat FS RT675 UP	
Receiver type	UP
Output / Gain (2cc) Filtered Earhook	134 / 78
Output / Gain Unfiltered Earhook	141 / 84




beat FS RT13 SP	
Receiver type	SP
Output / Gain (2cc) Filtered Earhook	133 / 75
Output / Gain Unfiltered Earhook	139 / 81



AQ beat FS R	
Receiver type	P
Output / Gain (2cc) Earhook	130 / 68
Output / Gain Slim Tube	129 / 69



beat FS RT312 M	
Receiver type	M
Output / Gain (2cc) Earhook	128 / 63
Output / Gain Slim Tube	122 / 58



--- Slim Tube



HANSATON is the exclusive brand for you and your clients since 1957.

We focus on personal relationships and your success.
Our German heritage, combined with the innovative power
of Sonova, enables us to create beautifully designed
products, based on cutting-edge, proven technology.



hansaton.com/pro



We drive the transition to a resource-efficient and
low-carbon future of our operations throughout the
lifecycle of our products and services.

sonova.com/sustainability

